




view **N**ow

BIOTEC Laboratories is pleased to announce the launch of
View Now Syphilis Rapid Devices

- Simple, rapid tests, easy to interpret
- Minimum equipment/training required
- Accurate results available in 10 minutes
- High sensitivity and specificity
- Detect IgM and IgG antibodies
- Storage at 2-30°C



www.biotec.com

BIOTEC Laboratories Ltd. 32 Anson Road, Martlesham Heath, Ipswich, Suffolk, IP5 3RG, UK.
Telephone: +44 (0) 1473 612158, Fax: +44 (0) 1473 611476
Email: sales@biotec.com, Website: www.biotec.com

view NOW Syphilis Rapid Devices

Possible applications include:

- Screening antenatal patients in remote clinics unable to perform rapid plasma reagin (RPR) tests.
- Screening patients that would otherwise be lost to follow-up.
- As a screening or diagnostic test in populations with a low rate of past, treated syphilis infection.
- Use as an on-site outbreak investigation tool.
- Field use in mass testing and treatment protocols (e.g. inmates, military recruits, returning migrant workers).
- Use in refugee settings and other complex emergency situations that interrupt access to routine testing.
- Enhancing clinical protocols for genital ulcer disease and to identify or rule out syphilis.
- Epidemiological surveillance in prevalence studies.
- Confirmation of non-treponemal test results in primary healthcare settings.
- Use in innovative settings such as antenatal/family planning clinics, pharmacies, private providers, pre-marital testing, and prison health clinics.

Ref.	Name	Sample Type	Tests per kit	Sensitivity*	Specificity*	Accuracy*
8/502	Syphilis Rapid Device	Serum/Plasma	40	>99.9%	99.7%	99.8%
8/504	Syphilis Rapid Device	Serum/Plasma/ Whole Blood	40	99.7%	99.6%	99.7%

*vs TPHA.

Also available from BIOTEC in the View Now Range:

hCG Rapid Strips in a tube	25 strips €€	Malaria Pf Rapid Device	25 test €€
hCG Rapid Strips in single pouches	50 strips €€	Malaria Pf/Pv Rapid Device	25 test €€
HIV Rapid Devices	40 test	Malaria Pan/Pf Rapid Device	25 test €€
HBsAg Rapid Devices	40 test	Malaria Pan/Pf/Pv Rapid Device	25 test €€
HCV Ab Rapid Devices	40 test		

Rev. 1 (2006-09)

www.biotec.com

BIOTEC Laboratories Ltd. 32 Anson Road, Martlesham Heath, Ipswich, Suffolk, IP5 3RG, UK.
 Telephone: +44 (0) 1473 612158, Fax: +44 (0) 1473 611476
 Email: sales@biotec.com, Website: www.biotec.com